



## SINGLE DISC POLISHING MACHINE PL-1



### **OVERVIEW :**

- For secondary surface preparation of metallurgical specimens in the laboratory, it is an effective & heavy duty machine with variable speed for polishing specimens for microscopic analysis.
- Both polish paper & polishing cloth can be clamped on rotating disc, provision for water flow & control.

### **FEATURES :**

- No. of Disc – One
- Disc diameter – 20 cms
- Power supply – Single Phase, 220 VAC, 50 Hz.
- Drive Motor – 0.5 HP with Variable frequency drive.
- Speed – Continuously variable upto 3000 RPM Maximum.
- Dimensions - L51XW35XH36 cms.

## VERTEX ENGINEERS & ASSOCIATES

R.H- 116, Sector 2B, Near Datta Meghe Engineering College, Airoli, Navi Mumbai - 400 708

Tel : +91-22-27791675 / 65169832 Email : [info@vertexmachines.com](mailto:info@vertexmachines.com), Visit us at : [www.vertexmachines.com](http://www.vertexmachines.com)