



VERTEX

AUTOMATIC SPECIMEN MOUNTING PRESS



The mounting press is equipped with a digit temperature controller to establish real-time temperature meanwhile, a timer is equipped to realize sample making Automatically increasing working efficiency. For each specimen-mount making 3-4 minutes is enough and a mirror-like specimen can be acquired, High speed mould heater with built-in water cooling coils. Interchangeable mould. Suitable for high volume application, relieves operator fatigue



SPECIFICATIONS:

- LCD 7" Touch Screen HMI & PLC Based Control Panel
- Heater Capacity : 1500 watt
- Automatic Water cooling systems
- Heating & Cooling Time : Adjustable 1-99 minutes
- Temp Range : 50-400 degree Adjustable
- Pressure : 50 -100 PSI Adjustable

- Inter Changeable Mould Size Dia-25, 30, 40, 50 mm
- Aluminum Casting Heater with Water Cooling Provision
- Alarm System after Completing the Cycle
- Pneumatic Operation System Compressed air required
- Quick ram return to speed up reloading
- Dovetail slot for easy changeability of Mould assemblies
- Digital temperature indicator and digital timer with LCD display
- Standard mould size 31.75 Or 38 mm dia provided with the machine
- Quick heater for speedy operation
- Pressure gauge, FRC
- Power Supply : 230 C/50 Hz
- Machine Size: 400X650X450 mm



<https://www.vertexmachines.com/>



+91 9323405785



info@vertexmachines.com

Actual product may differ in color and appearance due to constant development.