



VERTEX

DENSITY BALANCE



Density Determination Kit is available in 2 types one to be mounted on the pan and the other type is with attachment for weigh below the balance. Density Determination Kit consists density weighing pan assembly, stainless steel base for beaker, glass beaker, optional sinker for testing materials having density less than water and software (program) for density determination.

- Facility for programming density of liquid (other than water) and temperature of water.
- In case of multiple samples, storage of all the sample weights in air first and then in water, viewing of density of individual samples one by one on the display.
- Determination of air release value for oil samples (Optional).

SPECIFICATION:

MODELS	DM-120	DM-220	DM-320	DM-410	DM-500
CAPACITY	120g	220g	320g	410g	500g
READABILITY	0.001g	0.001g	0.001g	0.001g	0.001g
REPEATABILITY (+/-)	0.001g	0.001g	0.001g	0.001g	0.002g
LINEARITY (+/-)	0.001g	0.002g	0.002g	0.002g	0.002g
PAN SIZE	100mm	100mm	100mm	100mm	100mm