



VERTEX

PROFILE PROJECTOR



SPECIFICATION:

Model	VM200
Screen Diameter	200 mm
Projector Screen	Rotatable 360° with 2 Minute vernier resolution
Magnification (Standard)	10X
Magnification (Optional)	20X, 25X, 50X, 80X
Magnification accuracy	+/- 0.05% on contour & +/- 0.01% on surface
Work Stage	Medium duty 350X200 mm fitted with high precision linear bearing
Stage travel X-Y	150 X 100 mm
Focusing Range	100 mm
Resolution	0.01 mm / 0.005 mm / 0.001 mm
Measuring Range	With mechanical micrometer 0-25 mm
Illumination	12V- 100W Halogen lamp with intensity controller
Filter Colour	Green
Power Supply	230V AC 50 Hz
Standard Accessories	1 No. Projection Lamp, Operation Manual



<https://www.vertexmachines.com/>



+91 93234 05785



info@vertexmachines.com

Actual product may differ in color and appearance due to constant development.