



# VERTEX

## SCUFF TESTER



### OVERVIEW:

This equipment is intended to evaluate the Rub Proofness of prints on paper or board. It can also be used to measure / evaluate Color Transfer from printed / coated surface during rubbing. The application can be extended for measurement of abrasion resistance of some Plastic materials / Aluminum Foils.

### FEATURES:

- The highest degree of parallelism of disc samples.
- Accurate R.P.M. of rotating disks assured by High torque Step Synchronous motors
- Very Low Noise Level.
- Compact but Sturdy Design
- Very Low power consumption (Less than 100w).
- User - Friendly Electronic Controls.
- Non - Contact Type Sensor for Rub Counter makes it maintenance free.

### SPECIFICATION:

Specimen Size	2 inch and 4.5 inches Diameter
Weight on Sample	2psi & 4psi separate weights
Controls	Digital Preset Cycle counter
No. of Cycles	0 - 9999, 9999
Paint	Powder coated & Chromate finished
Input Power	230 Volts, 50/60Hz, single phase supply

**DEMO VIDEO: Please Click on it to open**

