



VERTEX

LEAK (VACUUM) TEST APPARATUS



This apparatus is an ideal instrument to check sealing integrity of Blister Packs; Strip Packs; Sachet and Pouch Packs etc. It is also used for checking leaks from Capped Bottles; Induction Sealed Containers; Sealed or Crimped Tubes etc. The old method of separate desiccator and huge vacuum pump is replaced with this new Leak Test Apparatus in GMP Stainless Steel body with phenomenal compactness.

SPECIFICATION:

Vacuum Range	Up to 550 mm of hg.
Desiccator Size	300 mm Dia.
Hold Time	5 Minutes

The instrument is equipped with long lasting vacuum pump, which gives higher level of vacuum in shortest possible time. Die pressed polycarbonate vacuum desiccator sustains vacuum for long time.

A vacuum gauge is provided to indicate vacuum level. A vacuum flow control cum isolation valve is provided to disconnect the vacuum source, thus avoiding the need of continuous running of vacuum pump. A digital timer is provided to set vacuum hold time as desired. The apparatus's amazing compactness saves lot of valuable space in process control laboratories.



<https://www.vertexmachines.com/>



+91 93234 05785



info@vertexmachines.com

Actual product may differ in color and appearance due to constant development.