



VERTEX

SINGLE DISC POLISHING MACHINE MODEL PL-1



OVERVIEW:

- For secondary surface preparation of metallurgical specimens in the laboratory, it is an effective & heavy duty machine with variable speed for polishing specimens for microscopic analysis.
- Both polish paper & polishing cloth can be clamped on rotating disc, provision for water flow & control.

FEATURES:

No. of Disc	One
Disc diameter	20 cms
Power supply	Single Phase, 220 VAC, 50 Hz.
Drive Motor	0.5 HP with Variable frequency drive.
Speed	Continuously variable upto 3000 RPM Maximum.
Dimensions	L51XW35XH36 cms.



<https://www.vertexmachines.com/>



+91 93234 05785



info@vertexmachines.com

Actual product may differ in color and appearance due to constant development.